

MOHAWK

COLLEGE OF APPLIED ARTS AND TECHNOLOGY

**CHEMICAL, ENVIRONMENTAL, AND BIOTECHNOLOGY
DEPARTMENT**

Spectrometry: Quantitative Determination of ASA by Absorbance of Visible Light

by Professor David Cash

September, 2008

**Mohawk College is the author and owner of these materials
(excluding copyright held by others) and all copyright and
intellectual property rights contained therein.**

**Use of these materials for teaching or other non-commercial
purposes is allowed.**

**Contact information for Mohawk College will be found on
the following page.**

This Experiment is a 3 hour Analytical Chemistry laboratory exercise. It is designed for students in a common second term course of a 2-year diploma program (Biotechnology, Environmental, or Health Technician).

**For Information or Assistance
Contact:**

**MOHAWK COLLEGE
CHEMICAL, ENVIRONMENTAL, AND
BIOTECHNOLOGY DEPARTMENT**

**Professor Cindy Mehlenbacher
cindy.mehlenbacher@mohawkcollege.ca
905-575-1212 ext. 3122**

**Bill Rolfe (Chief Technologist)
bill.rolfe@mohawkcollege.ca
905-575-2234**

Experiment 8

Spectrophotometry: Analysis of the ASA content of a Tablet by Use of the Beer-Lambert Law

OBJECTIVE

A spectrophotometric analysis will be performed. The comprehension and skills learned will be transferable to other laboratory and workplace situations.

- A set of ASA solution standards will be prepared; the photometric absorbances of the solution standards at **530 nm** wavelength will be used to construct a least squares linear calibration curve.
- The calibration curve will be used to estimate the ASA content of an ASA tablet unknown.

REFERENCE

Harris, Chapter 4, pages 79-85 and Chapter 18, pages 378-390.

INTRODUCTION

Spectrophotometry

This method of analysis, sometimes called **spectrometry**, refers to a method based on the measurement (**metry**) of light energy (**photo**) at a selected wavelength (**spectro**).

Harris (Chapter 18) gives an introduction to the theory and practice of spectrometry.

What follows is a practical introduction to using **molecular absorption spectrometry** in an analysis.

Molecular Absorption Spectrometry

Some molecular substances and some polyatomic ions in solution absorb ultraviolet, visible, or infrared light energy. The amount of the absorption of the light varies with the substance, with the selected wavelength of the light, and with the temperature.

Quantitative Spectrometry

Lambert (~ 1760) and **Beer** (~ 1850) determined that the transmission of light energy through a solution of an absorbing substance was dependent on the **path length** and the **concentration** of the absorbing substance by a logarithmic relationship, as shown:

Beer - Lambert Law

(Beer's Law):


$$A = abc$$

A = the Absorbance

a = the 'Constant'

b = the Path Length

c = the Concentration



P_0 = the Incident Intensity

P = the Emergent Intensity

Definitions and Terms:

$$\text{Transmittance} = T = \frac{P}{P_0}$$

$$\text{Percent T} = \% T = \frac{P}{P_0} \times 100 \%$$

$$\text{Absorbance} = A = -\log T = \log \frac{P_0}{P} = \log \left(\frac{100}{\% T} \right)$$

This relationship is now called the **Beer-Lambert Law** or **Beer's Law**.

The Calibration Curve Method

Most quantitative methods using Beer's Law are based on the premise that the curve of absorbance against concentration of the analyte is **linear** in the region of use. This must be proven for a new method, and should be verified for an analysis that is unfamiliar.



The following data is taken from **Skoog**¹.

The data represents the absorbance at some unspecified wavelength, of solutions of Glucose.

| Concentration (mmol / L) | 0.02 | 2.00 | 4.00 | 6.00 | 8.00 | 10.00 |
|--------------------------|-------|-------|-------|-------|-------|-------|
| Absorbance | 0.002 | 0.150 | 0.294 | 0.434 | 0.570 | 0.704 |

On the left below is an **X - Y Scatter Plot** of the data with **NO LINES**. To the right is the same plot, to which a **TRENDLINE - LINEAR** has been added. The line equation is of type:

$$y = mx + b$$



This Section Continues on the Next Page →

¹ Skoog, West, Crouch and Holler, Fundamentals of Analytical Chemistry, 8th edition, page 221, Question 8-16.

The Calibration Curve Method (Cont.)

The equation of the line and the R^2 value are displayed on the plot.

The closer the value of R^2 is to **1.0000**, the better the fit of the data is to a straight line.

For this particular set of data the **slope** and **intercept** value of the linear trendline are very small numbers. In such a case, when the equation is displayed and selected, it is necessary to format the numbers in the equation to be displayed to a chosen number of decimal places. In this case, four (4) decimal places has been chosen.

Using Beer's Law in Analysis

Beer's Law can be used in an analysis in many ways. Starting from the best possible method:

1. The best way is to prepare a calibration curve from a set of **independently prepared** calibration standard solutions. This means each calibration solution is prepared from a separate weighing operation. Solutions of unknowns will have their concentration determined by using the calibration curve equation.
2. The second best way is to prepare a calibration curve from a set of calibration standard solutions prepared from a **single** stock solution. Solutions of unknowns will have their concentration determined by using the calibration curve equation. This method will fail if the stock solution is not properly prepared or if the dilutions are faulty, but this may not be immediately apparent. This is a normal method for experienced analysts, since all mass values and dilutions are assumed to be reliable, and quality assurance protocols are usually in effect to detect any problems in methodology.
3. A slightly less reliable method is **standard addition**. A known amount of analyte is added into the unknown solution. The absorbance is measured with and without the addition. This can work well if you are certain that you are in a linear region of the calibration curve.
4. The least reliable method is by **proportion**. The absorbance of the unknown solution is related to the absorbance of a single calibration standard solution. This can work well if you are certain that you are in a linear region of the calibration curve. In this method, the following relationship is utilized:

$$\frac{\text{Concentration of Unknown}}{\text{Concentration of Standard}} = \frac{\text{Absorbance of Unknown}}{\text{Absorbance of Standard}}$$

Uncertainty in Using a Calibration Curve Method

Skoog (8th edition, pages 207-208) discusses the uncertainty of using a calibration curve method in detail. **Harris** covers most of the same points without being as explicit. The uncertainty of the slope and intercept values lead to increasing uncertainty at the outer ends of the calibration region.

When using a calibration curve with a linear trendline, always keep the unknowns at the centre of the linear calibration curve range.

Precision and Accuracy in Spectrophotometry

Harris discusses the precision and accuracy of using a spectrophotometer (page 385) and advises keeping the Absorbance of the solutions between the values **0.4** and **0.9**.

Skoog, 8th Edition gives a more detailed discussion (Figure 26-11, page 801); the accuracy and precision of every instrument is different. A **Spectronic 20**, the type of instrument you will probably use, is accurate and precise to within **2 %** over the range of Absorbance from about **0.2** to **0.9**. More expensive instruments can be accurate and precise to within 1 % or better over a range of Absorbance from 0.1 to 2.0 or greater.

Constraints on a Beer's Law Analysis

To use Beer's Law in analysis **successfully**, some **caution** is required.

- Keep the **temperature constant** (room temperature is usually used).
- Keep the **wavelength fixed**, and as narrow a range as possible.
This depends on the nature and quality of the instrument used.
- Keep the cell path length constant (use the **same cell** for all measurements).
- Keep solutions **free of dust or other solids**, which block or scatter the transmission of light.
- Keep the cell surface **clean, and free of dirt, grease, and scratches**.
- The absorbance curve may have a limited linear range of concentrations.
Keep to the linear range.
- Instrument error increases at very low absorbance and at very high absorbance.
This depends on the instrument. See **Harris**, Figure 18-8, page 385.
Keep to the middle range of absorbances (0.2 – 0.9) for all analysis solutions.
- Check for interference. If there is some other substance present that absorbs at the same wavelength, called a 'matrix' effect, this will increase error greatly.
Check for background absorbance at the analytical wavelength.

Acetylsalicylic Acid


The beneficial properties of **salicylic acid** (2-hydroxybenzoic acid) have probably been known to human beings for many tens of thousands of years. This compound, widely found in nature in the roots, bark, leaves and fruits of many plants and trees, is an **analgesic** (relieves pain), an **antipyretic** (reduces fever) and an **anti-inflammatory** (reduces swelling). It can be taken orally (by mouth), but is very acidic and irritates the stomach lining severely.

The compound **acetylsalicylic acid**, also known as ASA or Aspirin[®], was first synthesized in 1853. This compound breaks down rapidly in the body to form salicylic acid, giving all the same beneficial effects as salicylic acid, but it can be taken orally with far less irritation of the stomach.

Acetylsalicylic acid is a solid at room temperature, existing as a white powder or white crystalline needles². It has a melting point of 135 – 137 °C, and decomposes at 140 °C. It is not very soluble in water, about 0.33 g / 100 mL at 25 °C, but is much more soluble in alcohols and acetone. Both salicylic acid and ASA are toxic and can be fatal in excess. The LD₅₀³ for acetylsalicylic acid⁴ in rats is about 200 mg / kg. The LD₅₀ for salicylic acid⁵ in rats is about 900 mg / kg.


Synthesis of Acetylsalicylic Acid

The commercial synthesis of **acetylsalicylic acid** utilizes the reaction of **salicylic acid** with **acetic anhydride**, using acid catalysis as shown in the equation below. The byproduct **acetic acid** is also recovered and reused.



ASA as a Carboxylic Acid

Since ASA is a carboxylic acid, it acts as a weak acid and can be neutralized by an alkali to form a salt. The K_a value of ASA is 3.3×10^{-4} .



² Merck Index, 9th Edition, 1976, Merck and Co., Entry 874, page 114.


³ LD₅₀ (lethal dose 50 %): the amount of a toxic agent (as a poison, virus, or radiation) that is sufficient to kill 50 percent of a population of animals usually within a certain time.

⁴ <http://www2.siri.org/msds/f2/cfg/cfgqj.html> (retrieved 2006 06 08)

⁵ <https://fscimage.fishersci.com/msds/20315.htm> (retrieved 2006 06 08)

Stability of Acetylsalicylic Acid

Acetylsalicylic acid is not chemically inert. If exposed to moisture, over time it breaks down to form salicylic acid and acetic acid. A sealed bottle of ASA tablets that is kept past its expiry date will often have an odour of acetic acid (vinegar) when opened:



The breakdown of ASA is much more rapid in the presence of acid or base, and at higher temperature. Bottles of low-dose ASA usually contain a small canister of **silica gel** mixed with **activated carbon**. The silica gel absorbs water vapour, and the activated carbon absorbs acetic acid vapour.

Solubility Considerations

ASA is not very soluble in water. Also, solid crystals of ASA are very slow to dissolve at room temperature. The ASA in tablets for human consumption are in the form of a very fine powder. The powdered ASA dissolves much more rapidly than large crystals of solid.

Tablet Coatings and Binders

ASA tablets sold as over-the-counter medications may have two kinds of coatings. Most tablets have a coating that makes them easy to swallow, but disintegrates rapidly in water. Some tablets are enteric⁶ coated. The enteric coating does not dissolve rapidly in the acidic solution present in the stomach, but does dissolve or soften in the neutral or slightly basic solution present in the small intestine. An enteric coating is used to prevent exposure of the stomach or esophagus to ASA. Some people find ASA to be a very irritating and dangerous drug.

The enteric coating on an ASA tablet is usually a polymeric substance. Although it softens and allows the powder inside to escape, it does not fully dissolve.


Most tablets contain binder substances, which hold the powdered material into a pill shape. Starch is often used as a binder. The binders are usually insoluble in water.

Coating and binder substances can be removed by a filtration step, to prevent interference with the analytical measurements.

⁶ **Merriam-Webster Dictionary: enteric** - being a coating (as of an aspirin tablet) designed to pass through the stomach unaltered and disintegrate in the intestines.

Analysis of ASA by Spectrometry

The analysis of ASA by spectrometry is actually the analysis of **salicylic acid** by spectrometry. Acetylsalicylic acid (ASA) is rapidly decomposed to salicylic acid by boiling in sodium hydroxide solution. It is a convenient fiction in the analysis to calculate solution concentrations as if the ASA were still present




Formation of the Ferric Ion Complex of Salicylic Acid


Ferric ion is well known for its ability to form intensely coloured complexes with phenols. This reaction is characteristic of almost all phenols, including salicylic acid. This reaction is used as a test for the presence of a phenol.

Complex Ion formed by Salicylic Acid and Ferric Ion

The complex ion formed may not be exactly as shown here. The unlabeled positions around the ferric ion are probably occupied by the oxygen atoms of water molecules. Regardless, in water solution, in the presence of an excess of ferric ions, salicylic acid reacts **quantitatively** to form an intensely coloured substance. To the human eye, the solution appears to be **purple-violet** in colour.



A **maximum** of absorbance of visible light by the complex ion occurs at a wavelength of **530 nm**. This is the wavelength which should be selected for photometry.



An Example Data Set and an Example of a Calibration Curve for the Experiment


A **single** sample of pure ASA was weighed out using an analytical balance, heated in **1.0 M NaOH** solution, and diluted with distilled water into a **250 mL** volumetric flask to form a stock solution. Portions of this stock solution of ASA were measured out from a buret into a set of five **(5) 50 mL** volumetric flasks. This is **Method 2** on page 3.

The solutions were all diluted to **50 mL** total volume with a solution of **0.02 M** ferric ion in **0.1 M HCl** solution. This latter solution was used as a blank or background to zero the instrument. All of the solution standards were **violet – purple** to the human eye. The absorbance of each solution was measured in duplicate at **530 nm** wavelength on a **Spectronic 20** instrument, using the same **1 cm** path-length test tube cell for all measurements.

The **Table** below gives the **concentration** and the **mean absorbance** of each of the five calibration solution standards as though the ASA were still present, a convenient fiction. The concentrations are given in **mg / L** units, also called **ppm** or **parts per million** units. An Excel X – Y Scatter plot of this data, the best-fit least-squares linear trendline, the line equation, and the **R²** value are shown.

For the data set of this experiment, it will be necessary to choose to display the slope and intercept values of the trendline equation to **five (5)** places after the decimal point.

| ASA Conc. (mg / L) | Mean Absorbance |
|-----------------------|--------------------|
| 15.9 | 0.097 |
| 31.8 | 0.192 |
| 48.4 | 0.292 |
| 65.6 | 0.396 |
| 84.0 | 0.508 |



Dilution

Dilution is a technical term used in the laboratory when the volume of a solution is increased. To dilute means to add solvent to increase the volume of a solution, so that the concentration of the solution goes from a higher value to a lower value.

Dilution is used very often in the laboratory. Many reagents are purchased as highly concentrated solutions which must be reduced in concentration for daily use. In instrumental analysis, solutions of very low concentration are often required as standards of comparison. These are made by dilution of standard solutions. Dilution of reagent solutions used in excess or approximate quantity may be done approximately, using graduated cylinders and beakers. Standard solutions must be diluted as accurately and precisely as possible, using volumetric apparatus.

Simple Dilution

The most common dilution is the single step or simple dilution. The calculation of a simple dilution procedure is done using the dilution equation:

$$C_{\text{dilute}} \times V_{\text{dilute}} = C_{\text{concentrated}} \times V_{\text{concentrated}}$$

This looks like the equation of a titration, because both sides represent a quantity. That quantity is the amount of dissolved substance, which does not change on dilution. So there is no balanced equation or mol to mol ratio to worry about in a dilution. The dilution equation may be rewritten in terms of **ratios**:

$$\frac{C_{\text{concentrated}}}{C_{\text{dilute}}} = \frac{V_{\text{dilute}}}{V_{\text{concentrated}}}$$

It may also be written in terms of a **dilution factor**:

$$C_{\text{dilute}} = C_{\text{concentrated}} \times \left(\frac{V_{\text{concentrated}}}{V_{\text{dilute}}} \right) \leftarrow \text{Dilution Factor}$$


When the volume of a solution is increased by a certain factor, say **10 mL** to **50 mL**, this is referred to as a **1 to 5** dilution or a **five-fold** dilution. In this case, the new concentration (C_{dilute}) is **1 / 5** of the old concentration ($C_{\text{concentrated}}$).

Always use your chemical common sense to check the result of a dilution calculation. (Does it compute?) Notice that you can use any concentration unit in a dilution calculation and any volume unit, as long as you do not change the units during the calculation.

Serial Dilution

In a serial dilution, one dilution follows another. Serial dilution makes it easy to prepare solutions of very low concentration.

For example: A stock standard solution of zinc chloride contains **1100 ppm (mg / L)** zinc ion. An intermediate standard solution is prepared by diluting **10.00 mL** of the stock solution to **100.0 mL** in a volumetric flask. A working standard solution is prepared by diluting **5.00 mL** of the intermediate standard solution to **100.0 mL** volume. What is the concentration of each of the diluted solutions?



The calculation can be done in two steps, each a simple dilution. It can also be done in one step for the final solution by multiplying the dilution factors together and using the ratio equation:


$$C_{\text{dilute}} = C_{\text{concentrated}} \times \left(\frac{V_{\text{concentrated}}}{V_{\text{dilute}}} \right)$$

In the above dilution scheme, the overall dilution factor is $(1 / 10) \times (1 / 20) = (1 / 200)$.

Multiple Dilution

This is a procedure carried out when a range of standard comparison solutions is needed for an analysis. In the example above, a set of solutions might have been prepared by taking **1.00, 2.00, 3.00, 4.00** and **5.00 mL** of the intermediate standard solution and diluting each portion to **100 mL**.

The working standard solutions would be **1.10 ppm, 2.20 ppm, 3.30 ppm, 4.40 ppm** and **5.50 ppm** zinc respectively.



Sample Calculations

Preparation of a Stock Standard Solution of ASA

Example 1

A precisely weighed sample of **0.2018 g** of analytical-grade pure solid ASA was boiled in **1.0 M** NaOH solution, following the instructions of this experiment. The resulting solution was washed quantitatively into a **250 mL (0.2500 L)** volumetric flask, diluted to volume with distilled water, and mixed completely. The resulting solution was a stock standard solution of ASA.

Calculate the concentration of the ASA in the stock standard solution in **mg / L (ppm)** units. State the value to 4 significant figures.

Answer

Use the definition of concentration in mg / L (ppm):

$$\text{Concentration in mg / L (ppm)} = \frac{\text{Mass Solute (mg)}}{\text{Volume Solution (L)}}$$

$$\text{Concentration in mg / L (ppm)} = 0.2018 \text{ g} \times \frac{1000 \text{ mg}}{1 \text{ g}} \times \frac{1}{0.2500 \text{ L}}$$

$$\text{Concentration of ASA} = \underline{\underline{807.2 \text{ mg / L (ppm)}}}$$

Preparation of a Diluted Solution of Standard ASA

Example 2

A **4.00 mL** sample of the stock standard solution of ASA of **Example 1** was pipetted into a **50.00 mL** volumetric flask and diluted to volume with **0.02 M** ferric ion solution in **0.1 M** HCl solution. This is a diluted standard solution of ASA for a calibration curve.

Calculate the concentration of the ASA in the diluted standard solution in **mg / L (ppm)** units. State the value to 3 significant figures.

Answer

Use the dilution relationship:

$$C_{\text{dilute}} = C_{\text{concentrated}} \times \frac{V_{\text{concentrated}}}{V_{\text{dilute}}} = 807.2 \text{ mg / L (ppm)} \times \frac{4.00 \text{ mL}}{50.00 \text{ mL}}$$

$$\text{Concentration of Diluted ASA} = \underline{\underline{64.6 \text{ mg / L (ppm)}}}$$

The Sample Calculations Section Continues on the Next Page →

Sample Calculations (Cont.)

Use of a Calibration Curve

An ASA tablet, nominal ASA content **325 mg**, was boiled in **1.0 M NaOH** solution, following the instructions of this experiment. The resulting mixture was filtered and washed quantitatively into a **250 mL (0.2500 L)** volumetric flask and diluted to volume with distilled water.

The resulting solution is the **stock solution** of the tablet unknown.

A **2.00 mL** sample of the stock solution of the tablet unknown was pipetted into a **50.00 mL** volumetric flask and diluted to volume with **0.02 M** ferric ion solution in **0.1 M HCl** solution.

This is the **diluted solution** of the tablet unknown.

The absorbance of this solution was measured on the **same instrument** as for the calibration data and curve on page 8, using the **same test tube sample cell**, at the **same time**, without changing the **wavelength setting**.

Example 3

The mean of two absorbance measurements of the diluted solution of the tablet unknown at **530 nm** wavelength was **A = 0.311**.

Use the calibration curve equation on page 8 to calculate the ASA concentration of the diluted solution of the tablet unknown in **mg / L** units. State the value to 1 place after the decimal point.

Answer

The concentration of the ASA could be estimated graphically from the calibration curve, as shown to the right. However, it is better to use the equation of the best-fit least-squares line.

$$y = 0.00604x + 0.00040$$

$$y = \text{mean A value} = \underline{0.311}$$

$$x = \text{ASA Concentration in mg / L units}$$


Substitute and solve:

$$0.311 = (0.00604)(\text{ASA Conc.}) + 0.00040$$

$$\text{ASA Conc.} = \frac{(0.311 - 0.00040)}{0.00604} = 51.4238... = \underline{51.4} \text{ mg / L}$$

Note: The value of the slope must be known to **5 places** after the decimal point, or the resulting calculated value will not be sufficiently precise.

The Sample Calculations Section Continues on the Next Page →



Sample Calculations (Cont.)

Example 4

The diluted tablet unknown solution was determined in **Example 3** to have an ASA concentration of **51.4 mg / L**. Use the dilution relationship to calculate the ASA concentration of the stock solution of the tablet unknown. The dilution was from **2.00 mL** volume of the concentrated solution to **50.00 mL** volume of the diluted solution. State the value to the nearest **1 mg / L**.

Answer

Rearrange the simple dilution equation on page 9:

$$\frac{C_{\text{concentrated}}}{C_{\text{dilute}}} = \frac{V_{\text{dilute}}}{V_{\text{concentrated}}}$$

$$C_{\text{concentrated}} = C_{\text{dilute}} \times \frac{V_{\text{dilute}}}{V_{\text{concentrated}}}$$

$$C_{\text{concentrated}} = 51.4 \text{ mg / L} \times \frac{50.00 \text{ mL}}{2.00 \text{ mL}}$$

$$C_{\text{concentrated}} = \underline{1285 \text{ mg / L}}$$

Example 5

Use the total volume of the stock solution of the tablet unknown (**0.2500 L**) to calculate the amount of ASA present in the original solution, which is the amount of ASA in the tablet. State the value to the nearest **1 mg**.

Answer

$$\text{ASA Content (mg)} = 1285 \text{ mg / L} \times 0.2500 \text{ L} = 321.25 \text{ mg}$$

$$\text{ASA Content} = \underline{321 \text{ mg}}$$

The Sample Calculations Section Continues on the Next Page →

Sample Calculations (Cont.)

Example 6

A weighed sample of analytical-grade ASA and an ASA tablet unknown were treated according to the experiment procedure. The data in the **Table** to the right were obtained. Calculate the ASA concentration in the diluted unknown solution by the method of proportion. State the value to the nearest **0.1 mg / L**.

| ASA Conc. (mg / L) | Mean Absorbance |
|--------------------|-----------------|
| 48.7 | 0.295 |
| Diluted Unknown | 0.314 |

Answer

This is **Method 4** on page 3. Rearrange the relationship as shown, substitute and solve:

$$\frac{\text{Concentration of Unknown}}{\text{Concentration of Standard}} = \frac{\text{Absorbance of Unknown}}{\text{Absorbance of Standard}}$$

$$\text{Concentration of Unknown} = \frac{\text{Absorbance of Unknown}}{\text{Absorbance of Standard}} \times \text{Concentration of Standard}$$

$$\text{Concentration of Unknown} = \frac{0.314}{0.295} \times 48.7 \text{ mg / L}$$

$$\text{ASA Concentration in the Diluted Unknown Solution} = 51.8366... = \underline{51.8} \text{ mg / L}$$

| | | | | | |
|------|--|-----|--|------------|--|
| Name | | Day | | Start Time | |
|------|--|-----|--|------------|--|

PRE-LABORATORY PREPARATION

To be completed before the laboratory session.

To be submitted before beginning the experiment (20 points).

Questions: Answer in the space provided. **Show work.**
Your Mohawk College ID Number is **nnnnVWXYZ**.

Calibration Curve Exercise (7 points)

Suppose that the data in the **Table** below were obtained by the method of this experiment. Use your Mohawk ID Number to determine and enter the value of each datum⁷.

| ASA Conc. (mg / L) | ASA Conc. (mg / L) | Mean Absorbance | Mean Absorbance |
|-----------------------|-----------------------|--------------------|--------------------|
| 15.V | | 0.09Z | |
| 31.W | | 0.19Y | |
| 48.X | | 0.29X | |
| 65.Y | | 0.39W | |
| 84.Z | | 0.50V | |

Instructions:

- Use a software package (e.g.: Microsoft Excel[®]) to display an X - Y Scatter Plot of the data. **Choose the no-lines option.**
- Title** the plot, **including your name and the date**, and **label** both axes.
- Choose a **linear** least-squares best-fit **trendline** to the data.
- Display** the **trendline** on your plot and also display the **equation** and **R²** value of the trendline.
- Use the cursor to **select** the equation box.
Format as a **number** and choose the option **five (5)** places after the decimal point.
- Attach a print-out of your completed plot to the **PRE-LABORATORY PREPARATION**.

The PRE-LABORATORY PREPARATION Continues on the Next Page →

⁷ **Datum** (Merriam-Webster Dictionary) - something used as a basis for calculating or measuring. Plural: data or datums.

PRE-LABORATORY PREPARATION (Cont.)

Preparation of a Stock Standard Solution of ASA

A precisely weighed sample of **0.20XY g** of analytical-grade pure solid ASA was boiled in **1.0 M** NaOH solution, following the instructions of this experiment. The resulting solution was washed quantitatively into a **250 mL (0.2500 L)** volumetric flask, diluted to volume with distilled water, and mixed completely. The resulting solution was a stock standard solution of ASA.

- Q-1. Calculate the concentration of the ASA in the stock standard solution in **mg / L (ppm)** units. State the value to 4 significant figures. (2 points) Show work.
See **Example 1** on page 11.

$$0.20XY \text{ g} = \underline{\hspace{2cm}} \text{ g}$$

$$\text{ASA Concentration of the Standard Stock Solution (mg / L)} = \underline{\hspace{2cm}} \text{ mg / L}$$

Preparation of a Diluted Solution of Standard ASA

A **3.00 mL** sample of the stock standard solution of ASA of **Example 1** was pipetted into a **50.00 mL** volumetric flask and diluted to volume with **0.02 M** ferric ion solution in **0.1 M** HCl solution. This is a diluted standard solution of ASA for a calibration curve.

- Q-2. Calculate the concentration of the ASA in the diluted standard solution in **mg / L (ppm)** units. State the value to 3 significant figures. (2 points) Show work.
See **Example 2** on page 11.

$$\text{ASA Concentration of the Diluted Stock Solution (mg / L)} = \underline{\hspace{2cm}} \text{ mg / L}$$

The PRE-LABORATORY PREPARATION Continues on the Next Page →

| | | | | | |
|-------------|--|------------|--|-------------------|--|
| Name | | Day | | Start Time | |
|-------------|--|------------|--|-------------------|--|

PRE-LABORATORY PREPARATION (Cont.)

Use of a Calibration Curve

An ASA tablet, nominal ASA content **325 mg**, was boiled in **1.0 M** NaOH solution, following the instructions of this experiment. The resulting mixture was filtered and washed quantitatively into a **250 mL (0.2500 L)** volumetric flask and diluted to volume with distilled water to prepare a stock solution.

A **2.00 mL** sample of the stock solution of the tablet unknown was pipetted into a **50.00 mL** volumetric flask and diluted to volume with **0.02 M** ferric ion solution in **0.1 M** HCl solution.

The mean of two absorbance measurements of the diluted solution of the tablet unknown at **530 nm** wavelength was **A = 0.31Z**.

- Q-3. Use **your** calibration curve equation from the **Calibration Curve Exercise** to calculate the ASA concentration of the diluted solution of the tablet unknown in **mg / L** units. State the value to 1 place after the decimal point. (5 points) Show work. See **Example 3** on page 12.

A = 0.31Z = _____

ASA Conc. the Diluted Stock Solution of the Unknown (mg / L) = _____ mg / L

The PRE-LABORATORY PREPARATION Continues on the Next Page →

PRE-LABORATORY PREPARATION (Cont.)

- Q-4. Use the dilution relationship to calculate the ASA concentration of the stock solution of the tablet unknown. The dilution was from **2.00 mL** volume of the concentrated solution to **50.00 mL** volume of the diluted solution. (2 points) State the value to the nearest **1 mg / L**. Show work. See **Example 4** on page 13.

ASA Conc. of the Stock Solution of the Unknown (mg / L) = _____ mg / L

- Q-5. Use the total volume of the stock solution of the tablet unknown (**0.2500 L**) to calculate the total amount of ASA in the tablet. (2 points) State the value to the nearest **1 mg**. Show work. See **Example 5** on page 13.

Total ASA Content of the Tablet Unknown (mg) = _____ mg

PRE-LABORATORY PREPARATION **Total =** / **20**

PROCEDURE

Ensure that the fume hood fans are switched **ON** and are operating.

There are no special disposal instructions for this experiment. All solids and solutions may safely be disposed of by way of the municipal solid waste containers or the sinks. However, **1.0 M NaOH** solution is **moderately corrosive** and **hazardous**, and must be used with **caution** and **respect**.

A. Preparation of Glassware and Apparatus

Work with a Partner

You will be assigned to work with a partner on **Parts A to E** of the experiment. Record the name of your partner on page 33 in the **DATA TABLES AND REPORT** section. The data from three to five partnerships will be combined as a supergroup to produce a calibration curve for the analysis.

The following **clean** glassware and laboratory apparatus is required for the experiment:

- | | |
|--|---|
| For each pair of students: | <input type="checkbox"/> a weighing bottle and its lid* |
| <input type="checkbox"/> a spatula | <input type="checkbox"/> two 250 mL volumetric flasks and their stoppers |
| <input type="checkbox"/> a glass stirring rod | <input type="checkbox"/> two 50 mL volumetric flasks and their stoppers (extra equipment) |
| <input type="checkbox"/> a long stem funnel | <input type="checkbox"/> a 10 mL Mohr graduated pipet |
| <input type="checkbox"/> three small beakers | <input type="checkbox"/> a rubber pipet squeeze bulb |
| <input type="checkbox"/> two erlenmeyer flasks | |

* The instructor may instead direct you to weigh the ASA using a clean, dry weighing boat.

- A-1. Clean the glassware and apparatus if necessary with a 1 % solution of detergent in warm water. See **Cleaning and Drying of Glassware** on page **Error! Bookmark not defined.** Rinse the cleaned glassware and apparatus with tap water and then with distilled water. To avoid breakage, do not leave any glassware standing in an unstable position.
- A-2. Dry the spatula, and the weighing bottle and its lid (if using the weighing bottle) in the oven at 110 or 120 °C for 15 minutes. **Carefully** remove the spatula and the weighing bottle and its lid from the oven on to a heat proof pad, and allow them to cool to room temperature before using them.
- A-3. The instructor will set up one or more hot-plates. Using large beakers, heat **25 mL** of distilled water per student for the filtration step in **Part B**. The filtration of Part B may be slow. Complete Part B as soon as possible.

Work with your Partner. Part B and Part C may be completed in sequence or simultaneously.


B. ASA Unknown – Stock Solution Preparation

- B-1. Label your two erlenmeyer flasks, to avoid confusion or loss later in a crowd of similar flasks. Label one to contain an ASA tablet unknown for analysis, and the other to contain a sample of analytical-grade ASA.
- B-2. Record the brand of your assigned ASA tablet unknown sample in **Table B**. Record a description of the tablet and the nominal ASA content (**mg**). Weigh the tablet accurately. Record the mass value in mg units in **Table B**. Place the tablet into its labeled flask.

Caution: 1.0 M sodium hydroxide solution is a **caustic** material. It is **corrosive** and **hazardous** to skin, eyes, clothing and materials. Wearing gloves would be a sensible precaution. Wash any spills off your skin immediately with cold water. A small spill (up to **25 mL**) may be collected into a damp cloth and flushed into the sink with cold water. For a large spill, consult the instructor immediately.

- B-3. **In the fume hood area**, use the supplied **10 mL** graduated cylinder or a dispenser as instructed to add about **10 mL** of **1 M NaOH** solution to each of the flasks.


- B-4. Rotate (swirl) the flask gently to mix well. Place the flask on one of the hot-plates set **in the fume hood**. The heat control must be on a **low** setting. Heat the solution to boiling. Remove **immediately** from the hot plate (**Caution: hot**). Do not boil away all of the liquid.



The ASA tablet unknown solution may have a residue. Tablets contain many non-medicinal ingredients, some of which are insoluble in water or NaOH solution.

- B-5. Allow the solution in the flask to cool until it is easily held. You may run cold water over the **outside** of the flask if you wish to speed the cooling.
- B-6. Label your two **250 mL** volumetric flasks. One is for the ASA standard preparation and the other for the ASA tablet unknown preparation.

- B-7. **ASA Tablet Unknown:** Use a clean long stem funnel and a wash bottle to **filter** and transfer the solution **quantitatively** from the beaker, into its clean **250 mL** volumetric flask. Follow the instructions on page **Error! Bookmark not defined.** to fold the filter paper (use fast filter paper).



- B-8. Rinse the beaker, the solid residue, the filter paper, and the funnel with three (3) small portions of hot distilled water, about **5 mL** each, adding the wash water to the volumetric flask.

Note: The solution may be slightly cloudy due to the use of soluble starch as a binder in the tablet. This will not interfere with the analysis.

This Section of the PROCEDURE Continues on the Next Page →


B. ASA Unknown – Stock Solution Preparation (Cont.)

- B-9. Add distilled water to the flask to about one cm below the mark line. Fill the flask to the mark line using a dropper pipet.
- B-10. Stopper the flask with a clean stopper. Hold the stopper in place with one hand. Turn the flask over **slowly** at least **17 times** to ensure that the solution is completely uniform.

C. ASA Standard – Stock Solution Preparation

- C-1. If not already done, **carefully** remove the spatula and the weighing bottle and its lid (if required) from the oven on to a heat proof pad, and allow them to cool to room temperature before using them.


- C-2. The desired or target mass of the analytical-grade ASA is **0.20 ± 0.02 g**. Place **0.18 g – 0.22 g** of the ASA into your clean, dry container using a **top-loading balance**. Use the method of weighing by addition from **Experiment 1**. Do not record these preliminary weighings.



Top-Loading Balance

Do not fill or empty your container over any type of balance.

- C-3. Use the method of weighing by difference from **Experiment 1** to weigh **accurately (analytical balance)** the sample of analytical-grade ASA into its clean labeled erlenmeyer flask. There will be some solid left in the container. This does not matter.




Analytical Balance

- C-4. Enter the mass values to four places after the decimal point in **Table C** in the **DATA TABLES AND REPORT** section. Determine and enter the value for the mass of ASA placed in the flask.

- C-5. **In the fume hood area**, use the supplied **10 mL** graduated cylinder to add about **10 mL** of **1 M NaOH** solution to the flask.

- C-6. Rotate (swirl) the flask gently to mix well. Place the flask on one of the hot-plates set **in the fume hood**. The heat control must be on a **low** setting. Heat the solution to boiling. Remove **immediately** from the hot plate (**Caution: hot**). Do not boil away all of the liquid.



The ASA standard solution should be clear and colourless, with no solid residue.

This Section of the PROCEDURE Continues on the Next Page →

C. ASA Standard – Stock Solution Preparation (Cont.)

C-7. Allow the solution in the flask to cool until it is easily held. You may run cold water over the **outside** of the flask if you wish to speed the cooling.

C-8. Use a clean long stem funnel and a wash bottle to transfer the solution **quantitatively** from the flask, into its clean **250 mL** volumetric flask.


C-8. Rinse the flask and the funnel with several small portions of distilled water, adding the wash water to the volumetric flask.

C-9. Add distilled water to the flask to about one cm below the mark line. Fill the flask to the mark line using a dropper pipet.

C-10. Stopper the flask with a clean stopper. Hold the stopper in place with one hand. Turn the flask over **slowly** at least **17 times** to ensure that the solution is completely uniform.

C-11. Calculate and enter in **Table C** the ASA concentration of the stock standard solution in **mg / L** units. See **Example 1** on page 11.

C-12. When **Part B** and **Part C** are completed, have your **Table C** initialed by the instructor.



Part D and Part E may be completed in sequence or simultaneously.

The stock solutions of the ASA unknown and the standard will be diluted and placed in an excess of ferric ion solution to develop the purple colour intensity which will be measured for the analysis.

D. ASA Unknown – Diluted Solution Preparation

D-1. **Label** your two **50 mL** volumetric flasks. One is for the **diluted** ASA tablet unknown solution. The other is for the **diluted** ASA standard solution.

Pipetting the ASA Unknown Stock Solution

D-2. **Label** one clean small beaker to be used for the **ASA tablet unknown stock** solution. Into this beaker, pour about **20 mL** of the ASA tablet unknown stock solution, using the beaker volume markings.

D-3. Rinse the inside walls of the beaker with the ASA tablet unknown stock solution. Use this portion of the solution to rinse out the **10 mL** graduated Mohr pipet as well. Collect and discard the rinsing solutions into the sink. Repeat the rinsing and discard the solution again. On the third refill, take about **10 mL** to **20 mL** of the ASA tablet unknown stock solution into the beaker.


For the ASA tablet unknown the objective is to have about **52 mg / L** of ASA in the final diluted solution. This will put the unknown absorbance approximately at the mid-point of the calibration curve as desired. The amount of the unknown stock solution to be taken depends on the nominal content of your tablet unknown. **You will probably be using a 325 mg tablet unknown.**

Dilute your ASA tablet unknown stock solution according to the appropriate entry in the **Table** following. **Circle** your assigned instruction here and also in **Table D**.

| | | | |
|---|-------------|-------------|-------------|
| Nominal ASA Content of Tablet (mg) | 81 | 325 | 650 |
| Volume of ASA Tablet Unknown Stock Solution (mL) | 8.00 | 2.00 | 1.00 |
| Nominal Final Concentration of ASA (mg / L) | 52 | 52 | 52 |

This Section of the PROCEDURE Continues on the Next Text Page →

The Mohr Graduated Pipet




Using a Mohr Graduated Pipet

A Mohr graduated pipet is used in a very similar manner to a buret, rather than a transfer pipet. The graduated stem of the Mohr pipet makes it possible to deliver non-standard and variable measured volumes. This makes the Mohr pipet very flexible. But it is more difficult to control and use than a buret.

Two Kinds of Mohr Pipets

A complication to the use of a Mohr pipet is the fact that there are **TWO** different kinds of Mohr pipet.

1. One kind of Mohr pipet is graduated in exactly the same manner as a buret (far left). The volume between the final graduation marking and the tip is “**undefined**” and unmeasured.
2. The other kind of Mohr pipet is more like a transfer pipet (near left). In this pipet, the volume between the final graduation marking and the tip is part of the delivered volume.



How to Use the Buret Type Mohr Pipet (At Left)


For example, if you require a **4.0 mL** delivery:

Fill the pipet to the **0.0 mL** line, and drain to the **4.0 mL** line. Or,

Fill the pipet to the **1.0 mL** line, and drain to the **5.0 mL** line. Etc.

However, you **cannot** fill the pipet to the **6.0 mL** line, and empty the pipet. The volume beyond the **10.0 mL** mark is **undefined**.

The buret type of Mohr pipet is called a **Measuring Mohr Pipet**



How to Use the Transfer Pipet Type Mohr Pipet (At Left)

For example, if you require a **4.0 mL** delivery:

Fill the pipet to the **0.0 mL** line, and drain to the **4.0 mL** line. Or,


Fill the pipet to the **1.0 mL** line, and drain to the **5.0 mL** line. Etc.

In this case, you **can** fill the pipet to the **6.0 mL** line, and empty the pipet. This pipet is designed to be emptied as part of the measured volume, making it easier to use and more precise in some applications.

The transfer pipet type of Mohr pipet is called a **Serological Mohr Pipet**

D. ASA Unknown – Diluted Solution Preparation (Cont.)

- D-4. Check that your squeeze bulb is clean and dry inside. Transfer by **10 mL Mohr pipet** the required volume (**2.00 mL** for a **325 mg tablet**) of **ASA tablet unknown stock solution** from its beaker into the appropriately labelled clean **50 mL** volumetric flask. **Wipe off the tip of the pipet before the transfer.** If you are unsatisfied with your pipetting technique in the transfer, discard the sample, rinse the volumetric flask with distilled water, and repeat.



**Never transfer by pipet directly from a volumetric flask or a storage bottle.
Always use a beaker or some other intermediate vessel.**

Dilution with Ferric Ion in HCl Solution

- D-5. Label one of the small beakers to be used with the **0.02 M ferric ion** solution (**Caution: 0.1 M HCl** in this solution is mildly hazardous). Into this beaker, pour about **20 mL** of the ferric ion solution solution, using the beaker volume markings.
- D-6. Rinse the inside walls of the beaker with the ferric ion solution. Use this portion of the solution to rinse out the **dropper pipet** as well. Collect and discard the rinsing solutions into the sink. Repeat the rinsing and discard the solution again. Refill the beaker as necessary to complete the preparation of the two diluted solutions (**Part D** and **Part E**).
- D-7. Add **ferric ion in HCl solution** to the **50 mL** volumetric flask. Mix well. Fill to about one cm below the mark line. Fill the flask to the mark line using a dropper pipet.
- D-8. Stopper the flask with a clean, **dry** stopper. Hold the stopper in place with one hand. Turn the flask over **slowly** at least **17 times** to ensure that the solution is completely uniform.

**Your diluted solutions should be lightly to moderately purple or violet in colour.
If it is not, immediately consult the instructor.**

E. ASA Standard – Diluted Solution Preparation

E-1. Transcribe the concentration of the stock standard ASA solution in **mg / L** units from **Table C** to **Table E** in the **DATA TABLES AND REPORT** section.

Pipetting the ASA Stock Standard Solution

E-2. **Label** one clean small beaker to be used for the **ASA stock standard** solution. Into this beaker, pour about **20 mL** of the ASA standard solution, using the beaker volume markings.


E-3. Rinse the inside walls of the beaker with the ASA standard solution. Use this portion of the solution to rinse out your **10 mL** graduated Mohr pipet as well. Collect and discard the rinsing solutions into the sink. Repeat the rinsing and discard the solution again. On the third refill, take about **10 mL** to **20 mL** of the ASA standard solution into the beaker.

The instructor will assign your partnership to be one of **three** to **five** pairs working to produce a set of data for a calibration curve.

You will dilute your ASA standard solution as one of the entries in the **Table** following. **Circle** your assigned number and volume here and also in **Table E**.

| Diluted Calibration Standard No. | 1 | 2 | 3 | 4 | 5 |
|--------------------------------------|------|------|------|------|------|
| Volume of ASA Standard (mL) | 1.00 | 2.00 | 3.00 | 4.00 | 5.00 |
| Target Concentration of ASA (mg / L) | 16 | 32 | 48 | 64 | 80 |


E-4. **Check that your squeeze bulb is clean and dry inside.** Transfer by **10 mL Mohr pipet** your assigned volume of ASA standard solution from its beaker into the appropriately labelled clean **50 mL** volumetric flask. **Wipe off the tip of the pipet before the transfer.** If you are unsatisfied with your pipetting technique in the transfer, discard the sample, rinse the volumetric flask with distilled water, and repeat.




ASA Standard Stock Solution (250 mL)



Intermediate Vessel



10 mL Mohr Graduated Pipet



Diluted ASA Calibration Solution (50 mL)

Never transfer by pipet directly from a volumetric flask or a storage bottle. Always use a beaker or some other intermediate vessel.

This Section of the PROCEDURE Continues on the Next Page →

E. ASA Standard – Diluted Solution Preparation (Cont.)

Dilution with Ferric Ion in HCl Solution

- E-5. Use the small beaker from **Part D** containing the **0.02 M ferric ion** solution and the dropper pipet (**Caution: 0.1 M HCl** in this solution is mildly hazardous).
- E-6. Add **ferric ion in HCl solution** to the **50 mL** volumetric flask. Mix well. Fill to about one cm below the mark line. Fill the flask to the mark line using the dropper pipet.
- E-7. Stopper the flask with a clean, **dry** stopper. Hold the stopper in place with one hand. Turn the flask over **slowly** at least **17 times** to ensure that the solution is completely uniform.

Your diluted solution should be lightly to moderately purple or violet in colour. If it are not, immediately consult the instructor.

F. Percent Transmission / Absorbance Measurements

The instructor will assign your partnership and larger group to a specific spectrophotometer. There should be a large beaker near the spectrophotometer for rinse solutions and other discarded solutions.

The instructor will demonstrate the preparation and use of the spectrophotometer. Each type of instrument used by the department has printed instructions for its use available in the laboratory. In general, the following steps are required:

- Plug in, turn on and allow sufficient warm-up time (about 10 minutes).
- Set the required wavelength (**530 nm**).
- Adjust **no transmission = 0 %** transmission of light = **infinite** Absorbance value.
- Adjust **full transmission = 100 %** transmission of light = **zero** Absorbance value.
These two adjustments are performed manually on the **Spectronic 20**, but are automatic on microprocessor equipped instruments.
- Clean your sample cell(s).
- Reset **full transmission** with either **distilled water** or a **blank solution** in the cell.

F-1 You will be given a special sample test tube, or be asked to share a sample test tube with (an)other student pair(s). Clean the test tube. Rinse well with distilled water after cleaning. Be careful not to scratch the surface of the glass.
You must use the same tube for all of your measurements.

Notice that the test-tube sample cell has a vertical white marker line at the top. This marker line is used to ensure that the tube is always placed into the instrument in the same position.

F-2. When your turn on the spectrophotometer comes, check that the wavelength is set to **530 nm**, but **DO NOT TOUCH** the wavelength setting dial if there is one.

F-3. Fill your sample test tube **three times** with the blank solution, the **0.02 M** solution of ferric ions in **0.1 M** HCl solution. Empty the blank solution out twice into a discard beaker. After the third refill, dry off the outside of the test tube carefully with the special low-lint wipes provided.

The tube needs to be only about two-thirds full; the light beam passes through the lower part of the tube.

F-4. Place the sample test tube in the sample compartment of the instrument with the white marker line in line with the marker notch on the cell compartment.
Reset **full 100 % transmission (Zero Absorbance)** with the **blank solution** in the test tube.

This Section of the PROCEDURE Continues on the Next Page →

F. Percent Transmission / Absorbance Measurements (Cont.)

- F-5. Check and reset the 0 % transmission reading with the cell compartment empty. **Be sure that the compartment is closed to outside light for this operation.** These settings interact, so you should repeat both settings at least twice more.
- F-6. Empty the sample test tube into the discard beaker. Rinse the test tube with either one of your prepared solutions **twice** and then fill the sample test tube a **third time** to be about two-thirds full. Dry the outside. Insert the test tube correctly into the instrument.

When the reading is stable, read and record the value for that solution in **Table F (Trial 1)** in the **DATA TABLES AND REPORT** section:

| Spectronic 20 | Other Direct Reading Instrument |
|---|---|
| Read % Transmission Scale Value to Nearest 1 % | Read Absorbance Value Displayed (All Digits) |

Note: DO NOT attempt to read the Absorbance scale of the **Spectronic 20** directly. The scale is not linear and is very difficult to interpret correctly.

- F-7. Repeat these steps for your other prepared solution. Record the value in **Table F (Trial 1)**.
- F-8. Repeat the process for both solutions if time permits (**Trial 2**). Do not worry if the values differ on the second reading by up to $\pm 2\%$ Transmission (or ± 0.01 Absorbance). If they differ by more than this amount, consult the instructor.
- F-9. Finally, repeat the process for the blank solution to check that the instrument reads **100 %** Transmission (or **0.000** Absorbance) or very close to that value with the blank solution in the light beam. Do not worry if the value differs from these values by up to $\pm 2\%$ Transmission (or ± 0.01 Absorbance). If it differs by more than this amount, consult the instructor.
- F-10. If you are using a **Spectronic 20**, you must **CALCULATE** the Absorbance of the two solutions as follows:
- Determine the mean value of % Transmission (= % **T**) for each solution. State the value to the nearest 1 %;
 - Calculate the value of Absorbance using this equation:

$$\text{Absorbance} = A = -\log T = \log \frac{P_0}{P} = \log \left(\frac{100}{\% T} \right)$$


The PROCEDURE Continues on the Next Page →

G. Data Entry (Work with your Group)

The instructor will combine your partnership with two, three, or four other partnerships as a supergroup. The data from all of the supergroup will be combined to prepare a calibration curve for the analysis which all of the pairs will be able to use for the report calculations.

Each supergroup will require a **recording secretary** and a **steering committee** to assess the suitability of each point of the data set. If time permits, unsuitable data points may be re-measured.

- G-1 Enter the student names of each of the partnerships of the supergroup in the appropriate cells in **Table G** in the **DATA TABLES AND REPORT** section.
- G-2. Enter the **ASA mass** and **Absorbance value** of the diluted solution standard for each of the partnerships of the supergroup in the appropriate cells in **Table G**.
- G-3. Calculate and enter the values of $C_{\text{concentrated}}$ and C_{dilute} in each appropriate cell of **Table G**.
- G-4. Plot the approximate positions of the values of **Absorbance** versus C_{dilute} for the data set. Use the graph area below.
- G-5. Show your **Table G** and the approximately plotted data values to the instructor for initialing and discussion. **If your data set is to be useful, the plot should be relatively linear.**



The REPORT Section is on the Next Page →

REPORT Steps R-1 to R-12 are to be carried out before leaving the laboratory.

Table B

R-1. Determine and enter the mass value of your ASA tablet unknown in **Table B**.

Table C

R-2. Enter your weighing data for the analytical-grade ASA in **Table B**.

R-3. Calculate the concentration of the ASA in the stock standard solution ($C_{\text{concentrated}}$) in **mg / L (ppm)** units. State the value to 4 significant figures. Show work. See **Example 1** on page 11. Enter the value in **Table C**.

Table D

R-4. Circle the volume of $C_{\text{concentrated}}$ that you used in the Table.

Table E

R-5. Circle the volume of $C_{\text{concentrated}}$ that you used in the Table.

R-6. Calculate the concentration of the ASA in the diluted standard solution (C_{dilute}) in **mg / L (ppm)** units. State the value to 3 significant figures. Show work. See **Example 2** on page 11. Enter the value in **Table D**.

Table F

R-7. Enter the experimental values of either **% Transmission** or **Absorbance** (depending on the instrument used) measured for each of your two solutions in **Table F**.

R-8. **Spectronic 20:** Determine and enter the mean value of % Transmission for each of the two solutions in **Table F**, and calculate and enter the Absorbance value for each solution as instructed. **Other Instruments:** Determine and enter the mean value of Absorbance for each of the two solutions in **Table F**.

Table G

R-9. Enter the student names of the partnerships in your calibration curve supergroup in **Table G**.

R-10. Enter the experimental data values of **ASA mass** and **Absorbance** values of the diluted standard solutions for all partnerships in your calibration curve supergroup in **Table G**.

R-11. Enter the calculated values of experimental $C_{\text{concentrated}}$ and C_{dilute} for each standard prepared in your calibration curve group in **Table G**.

R-12. Verify that all of the experimental data points for the calibration set are acceptable by preparing a rough plot of the data. Re-calculate or if necessary re-measure any points which are not correctly situated, as time permits.

The REPORT Section Continues on the Next Page →

REPORT (Cont.)

Table G (Cont.)

The remaining steps are to be carried out after the laboratory period.

R-13. Produce a calibration curve for the analysis.

Plotting Instructions (may be done either individually or in pairs)

- Use a software package (e.g.: Microsoft Excel[®]) to display an X - Y Scatter Plot of the data. **Choose the no-lines option.**
- Title** the plot, **including your name(s) and the date**, and **label** both axes.
- Choose a **linear** least-squares best-fit **trendline** to the data.
- Display** the **trendline** on your plot and also display the **equation** and **R²** value of the trend-line.
- Use the cursor to **select** the equation box.
Format as a **number** and choose the option **five (5)** places after the decimal point.
- Attach a print-out of your completed plot to the **REPORT**.

Analysis Calculations

R-14. Use the experimental calibration curve equation to calculate the experimental ASA concentration (C_{dilute}) of the diluted solution of **YOUR** tablet unknown in **mg / L** units. State the value to 1 place after the decimal point. Show work. See **Example 3** on page 12. Enter the value in **Analysis Calculations**.

R-15. Use the dilution relationship to calculate the experimental ASA concentration ($C_{\text{concentrated}}$) of the stock solution of **YOUR** tablet unknown. For a 325 mg tablet, the dilution was from **2.00 mL** volume of the concentrated solution to **50.00 mL** volume of the diluted solution. State the value to the nearest **1 mg / L**. Show work. See **Example 4** on page 13. Enter the value in **Analysis Calculations**.

R-16 Use the total volume of the stock solution of the tablet unknown (**0.2500 L**) to calculate the experimental total amount of ASA in **YOUR** tablet. State the value to the nearest **1 mg**. Show work. See **Example 5** on page 13. Enter the value in **Analysis Calculations**.

Bonus Questions (10 Points)

See the **Bonus Questions** on page 39.

| | | | | | |
|------|--|-----|--|------------|--|
| Name | | Day | | Start Time | |
|------|--|-----|--|------------|--|

DATA TABLES AND REPORT

The following **Tables** are to be used for recording observations and measurements. Measured values are to be recorded in the heavily shaded cells **IN INK**. Leave this page open on your bench during the experiment period. Have it initialed by the instructor on completing each section.

The completed and initialed **DATA TABLES AND REPORT** Section must be handed in along with any additional pages you may submit as your report.

Partner's Name: _____

Table B: ASA Unknown – Stock Solution Preparation

Description of the ASA Tablet Unknown

| Brand Name and Description of the ASA Tablet Unknown | Tablet Mass (mg) | Nominal ASA Content (mg) |
|--|------------------|--------------------------|
| | | |

Table C: ASA Standard – Stock Solution Preparation

ASA Target Mass: 0.20 ± 0.02 g

Weighings of Analytical-Grade Solid ASA (In Ink to 4 Places after the Decimal Point)

| Container + Solid (g) | Container + Residue (g) | Mass of ASA (g) |
|-----------------------|-------------------------|-----------------|
| | | |

Instructor's Initials: _____ (5 points)

Calculate the experimental concentration of the ASA in **YOUR** stock standard solution in **mg / L (ppm)** units. State the value to 4 significant figures. (2 points) Show work.
See **Example 1** on page 11. **Volume of Flask Used:** 250 mL = 0.2500 L

Experimental ASA Conc. of **YOUR** Stock Standard Solution = _____ mg / L

The DATA TABLES AND REPORT Section Continues on the Next Page →

DATA TABLES AND REPORT (Cont.)

Table D: ASA Unknown – Diluted Solution Preparation

Circle Your Assigned Volume of $C_{\text{concentrated}}$ in the Table.

| | | | |
|--|------|------|------|
| Nominal ASA Content of Tablet (mg) | 81 | 325 | 650 |
| Volume of ASA Tablet Unknown Stock Solution (mL) | 8.00 | 2.00 | 1.00 |
| Nominal Final Concentration of ASA (mg / L) | 52 | 52 | 52 |

Table E: ASA Standard – Diluted Solution Preparation

Circle Your Assigned Solution Number and Volume of $C_{\text{concentrated}}$ in the Table.

| | | | | | |
|--------------------------------------|------|------|------|------|------|
| Diluted Calibration Standard No. | 1 | 2 | 3 | 4 | 5 |
| Volume of ASA Standard (mL) | 1.00 | 2.00 | 3.00 | 4.00 | 5.00 |
| Target Concentration of ASA (mg / L) | 16 | 32 | 48 | 64 | 80 |

Concentration of the Stock Solution Standard (From Table C) = _____ mg / L

Calculate the experimental concentration of the ASA in **YOUR** diluted standard solution in **mg / L (ppm)** units. State the value to 3 significant figures. (2 points) Show work.
See **Example 2** on page 11. **Volume of Flask Used: 50.00 mL**

Experimental ASA Concentration of **YOUR** Diluted Standard Solution = _____ mg / L

The DATA TABLES AND REPORT Section Continues on the Next Page →

| | | | | | |
|------|--|-----|--|------------|--|
| Name | | Day | | Start Time | |
|------|--|-----|--|------------|--|

DATA TABLES AND REPORT (Cont.)

Table F: Percent Transmission / Absorbance Measurements

If Using the Spectronic 20

Percent Transmission Values (In Ink) to the Nearest 1 %

| | Trial 1 | Trial 2 | Mean % T |
|--|---------|---------|----------|
| % Transmission Value of Diluted Standard | | | |
| % Transmission Value of Diluted Tablet Unknown | | | |

Calculated Absorbance Values to 3 Places

| | Mean % T | Calculated Absorbance |
|------------------------|----------|-----------------------|
| Diluted Standard | | |
| Diluted Tablet Unknown | | |

$$\text{Absorbance} = A = -\log T = \log \frac{P_0}{P} = \log \left(\frac{100}{\% T} \right)$$

Example: Mean % T = 35 %

$$\text{Absorbance} = \log \left(\frac{100}{35} \right) = \underline{0.456}$$

Show a Sample Calculation

If Using a Direct Reading Instrument

Absorbance Values (In Ink) to 3 Figures

| | Trial 1 | Trial 2 | Mean A |
|--|---------|---------|--------|
| Absorbance Value of Diluted Standard | | | |
| Absorbance Value of Diluted Tablet Unknown | | | |

Instructor's Initials on Completion: _____ (10 points)

The DATA TABLES AND REPORT Section Continues on the Next Page →

DATA TABLES AND REPORT (Cont.)

Table G: Data Entry

| | Partnership Names |
|------------|-------------------|
| Standard 1 | _____ and _____ |
| Standard 2 | _____ and _____ |
| Standard 3 | _____ and _____ |
| Standard 4 | _____ and _____ |
| Standard 5 | _____ and _____ |

Calculated Concentrations and Absorbance Values of the ASA Solution Standards

| | Standard 1 | Standard 2 | Standard 3 | Standard 4 | Standard 5 |
|---|----------------|----------------|----------------|----------------|----------------|
| ASA Mass (g) | | | | | |
| Calculated ASA C _{concentrated} (mg / L) | | | | | |
| ASA Target C _{dilute} (mg / L) | 16 | 32 | 48 | 64 | 80 |
| Dilution of ? mL to 50.00 mL | <u>1.00</u> mL | <u>2.00</u> mL | <u>3.00</u> mL | <u>4.00</u> mL | <u>5.00</u> mL |
| Calculated ASA C _{dilute} (mg / L) | | | | | |
| Mean Experimental Absorbance | | | | | |

Instructor's Initials on Completion: _____ (10 points)

The DATA TABLES AND REPORT Section Continues on the Next Page →

| Name | | Day | | Start Time | |
|------|--|-----|--|------------|--|
|------|--|-----|--|------------|--|

DATA TABLES AND REPORT (Cont.)

Produce a calibration curve for the analysis.

Plotting Instructions (may be done either individually or in pairs) (7 points)

- Use a software package (e.g.: Microsoft Excel[®]) to display an X - Y Scatter Plot of the data. **Choose the no-lines option.**
- Title** the plot, **including your name(s) and the date**, and **label** both axes.
- Choose a **linear** least-squares best-fit **trendline** to the data.
- Display** the **trendline** on your plot and also display the **equation** and **R²** value of the trendline.
- Use the cursor to **select** the equation box.
Format as a **number** and choose the option **five (5)** places after the decimal point.
- Attach a print-out of your completed plot to the **REPORT**.

Analysis Calculations

Use the **experimental** calibration curve equation of **YOUR** group to calculate the experimental ASA concentration (**C_{dilute}**) of the diluted stock solution of **YOUR** tablet unknown in **mg / L** units.

State the value to 1 place after the decimal point. (5 points) Show work.

See **Example 3** on page 12.

Experimental ASA Conc. (**C_{dilute}**) of the Diluted Stock Solution of **YOUR** Unknown (mg / L)
= _____ mg / L

The DATA TABLES AND REPORT Section Continues on the Next Page →

DATA TABLES AND REPORT (Cont.)

Analysis Calculations (Cont.)

Use the dilution relationship to calculate the experimental ASA concentration ($C_{\text{concentrated}}$) of the stock solution of **YOUR** tablet unknown. For a 325 mg tablet, the dilution was from **2.00 mL** volume of the concentrated solution to **50.00 mL** volume of the diluted solution.

State the value to the nearest **1 mg / L**. (2 points) Show work. See **Example 4** on page 13.

Experimental ASA Conc. ($C_{\text{concentrated}}$) of the Stock Solution of **YOUR** Unknown (mg / L)
= _____ mg / L

Use the total volume of the stock solution of the tablet unknown (**0.2500 L**) to calculate the experimental total amount of ASA in **YOUR** tablet. State the value to the nearest **1 mg**. (2 points)

Show work. See **Example 5** on page 13.

Experimental Total ASA Content of **YOUR** Tablet Unknown (mg) = _____ mg

The Bonus Questions and the Mark Sheet are on the Next Two Pages →

| Name | | Day | | Start Time | |
|------|--|-----|--|------------|--|
|------|--|-----|--|------------|--|

DATA TABLES AND REPORT (Cont.)

Bonus Questions (10 Points)

Attach a separate page of calculations.

1. Calculate the ASA concentration (C_{dilute}) in **YOUR** diluted tablet unknown solution **by the method of proportion**. Use only your own data for **YOUR** diluted standard and **YOUR** diluted tablet unknown solution. State the value to the nearest **0.1 mg / L**. (6 points)
Show work. See **Example 6** on page 14.
2. Use the dilution relationship to calculate the ASA concentration ($C_{\text{concentrated}}$) of the stock solution of **YOUR** tablet unknown **using the value of C_{dilute} calculated as the answer to Bonus Question 1**. For a 325 mg tablet, the dilution was from **2.00 mL** volume of the concentrated solution to **50.00 mL** volume of the diluted solution.
State the value to the nearest **1 mg / L**. (2 points)
Show work. See **Example 4** on page 13.
3. Use the total volume of the stock solution of the tablet unknown (**0.2500 L**) to calculate the total amount of ASA in **YOUR** tablet unknown **using the value of $C_{\text{concentrated}}$ calculated as the answer to Bonus Question 2**. State the value to the nearest **1 mg**. (2 points) Show work. See **Example 5** on page 13.

The Mark Sheet is on the Next Page →

DATA TABLES AND REPORT (Cont.)

Mark Sheet for Experiment 8

Total = 100 points

| Category | Points |
|--|--------------------------|
| Pre-Laboratory Preparation | / 20 |
| Attendance: Punctuality, Diligence, Clean-Up. | / 25 |
| Data Acquisition and Recording In Ink /5 Proper Number of Digits /5 Results Initialed ASA Mass /5 Absorbances /10 Data Entry /10 | / 35 |
| Calibration Curve, Analytical Calculations | / 20 |
| Absorbance Value for the ASA Standard Solution Result Correct Within: $\pm 0.010\% = 20$ $\pm 0.019 = 16$ $\pm 0.036 = 12$ $\pm 0.069 = 8$ $\pm 0.130 = 4$ | / 20 |
| Bonus Questions Question 1 / 6 Question 2 / 2 Question 3 / 2 | / 10 |
| Comments: | Total = /120 |